

Catalase

Tests for the presence of the enzyme catalase.

Intended Use:

- Differentiates *Streptococcus* (-) from *Micrococcus* (+)
- Differentiates *Staphylococcus* (V+) and *Bacillus* (+) from *Clostridium* (-)

Principle:

Hydrogen peroxide (H_2O_2) is the end product of aerobic breakdown of sugars. Since it is toxic to bacterial cells, most aerobic bacteria produce catalase or peroxidase to protect themselves.

Streptococcus, *Enterococcus*, and *Lactobacillus* are exceptions. Since they do not use the cytochrome c pathway, they do not produce H_2O_2 and lack catalase.

Test Procedure:

1. Transfer a well isolated colony to a clean glass slide and add 1 drop of 3% H_2O_2 .
 - Do not reverse the order and do not mix.
2. Observe for immediate bubble formation.
3. Use 15% H_2O_2 for the detection of catalase in anaerobes.

Results:

- The formation of bubbles is considered a positive result.

Limitations:

- Do not take your colony from a blood agar plate. The catalase present in the erythrocytes will give a false positive result.
- H_2O_2 is unstable. You can do a quality control test of the H_2O_2 reagent by placing a drop on a blood agar plate. Vigorous bubbling should result.