

## Triple Sugar Iron (TSI) Agar Slant

### Principle:

To determine the ability of an organism to attack a specific carbohydrate incorporated in a basal growth medium, with or without the production of gas, along with the determination of possible hydrogen sulfide (H<sub>2</sub>S) production.

### Purpose:

Hydrogen sulfide production and/or fermentation patterns are generally characteristic for specific bacterial groups, genera, or species, especially among the Enterobacteriaceae.

### Test Procedure:

1. With an inoculating needle, stab the TSI slant about 2/3 of the way into the butt, withdraw the needle and streak the slant.
2. Incubate aerobically at 35-37°C.
3. Examine after 18-24 hours of incubation.

### Interpretations:

#### Carbohydrate utilization...

	<b>Fermentation of only glucose</b>	<b>Fermentation of glucose &amp; lactose</b>	<b>Neither glucose nor lactose fermented</b>
<b>Slant</b>	Alkaline reaction red color	Acid reaction yellow color	Growth only No color change
<b>Butt</b>	Acid reaction yellow color	Acid reaction yellow color	Growth only No color change

#### Gas Production...

- Aerogenic Gas production (CO<sub>2</sub>) evident by one of the following...
  - A single gas bubble or bubbles in the medium
  - Splitting of the medium
  - Complete displacement of the medium from the bottom of the tube leaving a clear area
  - Slight indentation of the medium from the side of the tube
- Anaerogenic - No gas production

#### H<sub>2</sub>S production:

The presence of black precipitate is evident by...

- A black color spread throughout the entire butt masking the acidity; may even be slight evidence on the slant
- A black ring near the top of the butt area
- A black precipitate scattered throughout the butt but not entirely

